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CLINICAL CASES: HEMI-ARTHROPLASTY FOR FRACTURE

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INSTALLATION AND APPROACH (CASE 1)

Insertion in a semi-seated position, with image intensifier facing towards the anaesthetists who can "arrive" during the operation and carry out the required monitoring (**fig. 1**).

The patient is placed on the edge of the table to allow a retro-pulsion of the limb to be operated on, without the head being in hyperextension: light flexion is preferable (**fig. 2**) position without causing difficulties in the procedure.

Preoperatively the direction of the image intensifier and the release of the articulating joint should be checked (**fig. 3**).

DELTO PECTORAL APPROACH

Cutaneous incision of coracoid in the axillary hollow without touching this latter (**fig. 4**). Unlike the supero external route, no muscle is cut and one enters "inside the fracture" above the largest vein. Continue on the side where there is the least collateral (**fig. 5**) or one can also immediately ligate the vein... which does not tolerate the placement of the retractor very well for the duration of the operation.

Then pass across the clavipectoral fascia (star) on the outside edge of the coracobrachialis (arrow) (**fig. 6**). Do not proceed into the medial, into the body of the coracobrachialis as this could create iatrogenic lesions on the musculocutaneous nerve.

Then one must find and isolate the head and tuberosities (**fig. 7**); Place markers on the length of the biceps which will not be cut straight away as it allows the symmetry between the tuberosities to be defined. The ridge is extracted and the tuberosities isolated.

Fig. 1

Fig. 2


Fig. 3


Fig. 4


Fig. 5


Fig. 6


Fig. 7

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3D reconstruction of the 4 part fracture

CASE 2

Female 64 yo, presenting with a CT4 fracture (**fig. 1**). The use of the humelock fracture prosthesis with the smart looped threads, and the OMS system with auto-transplant is followed by a useful (figure 2) and long-lasting (figure 3) functional result, because of a 1 month-long elbow to body immobilisation. At a review at 2.5 years, function is stable.

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Fig. 1: CT2 Fracture... in this age group it is often risky to suggest an osteosynthesis as a solution.

CASE 3

Female 73 yo, presenting with a CT2 fracture (**fig. 1**). Useful functional result after 1 year (**fig. 2**) with tuberosities consolidation (**fig. 3**).

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Fig. 3: X-ray view of tuberosity consolidation after 1 year... here one can see the how well the consolidation has taken, using the autograft and despite the less than optimal age group for osteogenesis.

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Fig. 1: *In case of 4 parts fracture between 60 and 70 yo, the implantation of hemiarthroplasty is a high demanding procedure.*

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Fig. 2: *Osteosuture of tuberosities and peroperative X ray control are mandatory in order to achieve bone union at the right place.*